



# Silver Staining of SDS-polyacrylamide Gel

Jennifer M. Kavran, Daniel J. Leahy<sup>1</sup>

Department of Biophysics and Biophysical Chemistry, Johns Hopkins University School of Medicine, Baltimore, MD, USA

<sup>1</sup>Corresponding author: e-mail address: dleahy@jhmi.edu

## Contents

1. Theory	170
2. Equipment	170
3. Materials	170
3.1 Solutions & buffers	170
4. Protocol	172
4.1 Duration	172
4.2 Preparation	172
4.3 Tip	172
4.4 Tip	172
4.5 Tip	172
4.6 Tip	172
5. Step 1 Fix the Gel	172
5.1 Overview	172
5.2 Duration	172
5.3 Tip	173
6. Step 2 Stain the Gel	173
6.1 Overview	173
6.2 Duration	173
7. Step 3 Preserve the Gel	175
7.1 Overview	175
7.2 Duration	175
References	176

## Abstract

To detect nanogram quantities of protein and nucleic acids on SDS-PAGE gels.

## 1. THEORY

Silver staining is more sensitive than both Coomassie blue (see [Coomassie Blue Staining](#)) and ethidium bromide staining (see [Agarose Gel Electrophoresis](#)). The chemistry of silver staining relies on the reduction of silver from an ionic state to a metallic state. Macromolecular bands are stained silver on a clear background. In practice, silver staining allows you to detect smaller amounts of protein and to easily check the purity of your sample.

## 2. EQUIPMENT

- Platform rotator
- Plastic gel staining box
- Gloves
- Gel documentation system
- Gel dryer or gel drying rack

## 3. MATERIALS

- Nitric acid ( $\text{HNO}_3$ )
- Ethanol, 95%
- Glacial acetic acid
- Glycerol
- Methanol ( $\text{MeOH}$ )
- Formaldehyde, 37% (Sigma)
- Sodium thiosulfate ( $\text{Na}_2\text{S}_2\text{O}_3$ )
- Sodium carbonate ( $\text{Na}_2\text{CO}_3$ , anhydrous, ACS grade)
- Silver Nitrate ( $\text{AgNO}_3$ )

### 3.1. Solutions & buffers

#### Step 1 Fixing solution

Component	Final concentration	Stock	Amount
Methanol	50%	100%	50 ml
Acetic acid	12%	100%	12 ml
Formaldehyde	0.0185%	37%	50 $\mu\text{l}$

Add water to 100 ml. Note: add formaldehyde immediately before use

Wash buffer

Mix 250 ml methanol with 250 ml water to give a final concentration of 50%

Step 2 Sodium thiosulfate solution

Dissolve 50 mg sodium thiosulfate,  $\text{Na}_2\text{S}_2\text{O}_3$ , in 250 ml water to give a concentration of 1.3 mM. Make fresh

Silver solution

Component	Final concentration	Stock	Amount
$\text{AgNO}_3$	11.8 mM	N/A	0.2 g
Formaldehyde	0.028%	37%	75 $\mu\text{l}$

Add water to 100 ml purified water. Make fresh

Developing solution

Component	Final concentration	Stock	Amount
$\text{Na}_2\text{CO}_3$ (anhydrous)	6% (w/v)	N/A	12 g
Sodium thiosulfate solution	0.05 mM	1.3 mM	4 ml
Formaldehyde	0.037%	37%	100 $\mu\text{l}$

Add water to 200 ml. Make fresh

Stop solution

Component	Final concentration	Stock	Amount
Methanol	50%	100%	50 ml
Acetic acid	12%	100%	12 ml

Add water to 100 ml

Step 3 Wash buffer 2

Dilute 30 ml methanol in 70 ml water to give a final concentration of 30%

Storage solution

Component	Final concentration	Stock	Amount
Ethanol	20%	100%	20 ml
Glycerol	10%	100%	10 ml

Add water to 100 ml



## 4. PROTOCOL

### 4.1. Duration

Preparation	1 h
Protocol	4 h

### 4.2. Preparation

Run SDS-PAGE mini-gel (see [One-dimensional SDS-Polyacrylamide Gel Electrophoresis \(1D SDS-PAGE\)](#)) and place in a clean gel staining box.

### 4.3. Tip

*Gel boxes should be cleaned with 50% nitric acid and then rinsed thoroughly with water.*

### 4.4. Tip

*Handle gels delicately. Wear gloves to avoid transferring fingerprints and avoid manipulating it with metal tools such as tweezers.*

### 4.5. Tip

*Staining will be uneven if the gel is not completely submerged in liquid. Increase buffer volumes as appropriate.*

### 4.6. Tip

*Times and volumes are appropriate for mini-gels. Increase both for thicker gels.*

See [Fig. 14.1](#) for the flowchart of the complete protocol.



## 5. STEP 1 FIX THE GEL


### 5.1. Overview

The proteins are immobilized and the SDS-PAGE running buffer is removed.

### 5.2. Duration

2.5 h

**1.1** Add 50 ml water to the gel in small container, swirl gently for 5 s, and decant water.



**Figure 14.1** Flowchart of complete protocol, including preparation.

- 1.2 Repeat Step 1.1 two to four more times.
- 1.3 Add 50 ml of Fixing Solution to the gel. Rotate the gel for at least 1 h. Carefully decant Fixing Solution.
- 1.4 Add 50 ml of Wash Buffer. Rotate for 20 min. Decant buffer.
- 1.5 Repeat Step 1.4 two more times.
- 1.6 Add 50 ml water to the gel in small container, swirl gently for 5 s, and decant water.
- 1.7 Repeat Step 1.6 two to four more times.

### 5.3. Tip

*The gel can be left overnight in Fixing Solution.*

See [Fig. 14.2](#) for the flowchart of Step 1.



## 6. STEP 2 STAIN THE GEL


### 6.1. Overview

The gel is incubated first in a silver solution and then in the developing solution, causing the silver to precipitate onto the protein.

### 6.2. Duration


1 h

- 2.1 Add 5 ml of sodium thiosulfate solution and incubate rotating for 1 min.



**Figure 14.2** Flowchart of Step 1.

- 2.2** Add 50 ml water to the gel in small container, swirl gently for 5 s, and decant water.
  - 2.3** Repeat Step 2.2 two to four more times.
  - 2.4** Add 50 ml of Silver Solution. Incubate rotating for 20 min. Decant solution.
  - 2.5** Add 50 ml water to the gel in small container, swirl gently for 5 s, and decant water.
  - 2.6** Repeat Step 2.5 two to four more times.
  - 2.7** Add 50 ml of Developing Solution. Quickly swirl in the gel box and decant.
  - 2.8** Add 50 ml of Developing Solution. Incubate rotating until desired staining intensity is achieved.
  - 2.9** Add 50 ml water to the gel in small container, swirl gently for 5 s, and decant water.
  - 2.10** Repeat Step 2.9 two to four more times.
  - 2.11** Add 50 ml of Stop Solution.
  - 2.12** Photograph the gel using a gel documentation system.
- See [Fig. 14.3](#) for the flowchart of Step 2.



**Figure 14.3** Flowchart of Step 2.

## 7. STEP 3 PRESERVE THE GEL


### 7.1. Overview

The gel is washed into a buffer for long-term storage.

### 7.2. Duration

45 min

**3.1** Decant the Stop Solution.



**Figure 14.4** Flowchart of Step 3.

**3.2** Add 50 ml of Wash Buffer (50% methanol). Incubate rotating for 20 min at 4 °C. Decant buffer.

**3.3** Add 50 ml of Wash Buffer 2 (30% methanol). Incubate rotating for 20 min at 4 °C. Decant buffer.

**3.4** Add 50 ml of Storage Solution. Incubate rotating for 20 min at 4 °C.

**3.5** The gel can now be dried using your preferred method.

See [Fig. 14.4](#) for the flowchart of Step 3.

## REFERENCES

Referenced Protocols in Methods Navigator

[Coomassie Blue Staining.](#)

[Agarose Gel Electrophoresis.](#)

[One-dimensional SDS-Polyacrylamide Gel Electrophoresis \(1D SDS-PAGE\).](#)